



- Insulated stainless steel cabinet
- Frost resistant down to (-40 °C)
- Additional frontal heating element (-20 °C to -40 °C)
- Fan for ventilation
- 2 x LED lights (*meaning defibrillator is visible in the dark*)
- Designed to be easy to open, even under conditions of heavy snow and frost.
- Future-proof! (*Software can be upgraded*)

FUNCTIONAL DESCRIPTION:

- If the door of the cabinet is opened... The SIREN immediately sounds.
(if the door of the cabinet is closed once again within 3 minutes the siren will stop automatically)
- Once the door is closed... The cabinet is once again ready for operation.
-

TECHNICAL DESCRIPTION:

- Door: The door is fitted with thermopane glass and marked with international AED labels/symbols.
 Light: 2 x LED lights in the cabinet, such that the defibrillator is visible in the dark.
 Siren: Siren in the cabinet. (100 dBA / 1m)
 Fan: Ventilation fan, starts automatically when temperature exceeds +40 °C.
 Frost-proof cabinet: Thermostatically-controlled heating element. Frost resistant down to -40 °C
 Voltage: Must be connected to a 230V (116Wh) electrical supply.
 Outdoor cable (3 x 1.5 mm²)
 Cat II, IE IEC 61010-1
- Overvoltage category: Cat II, IE IEC 61010-1
 Protection rating: IP43.
 Ambient temperature: -40 / +40°C
 Relative humidity: 0 – 95 %
 Environmental pollution degree: 2
 Altitude: Max. 2000 masl. (*meters above sea level*)
 Design norm: IEC61010 – with deviations in insulation fire class and service switch.
 Insulation – fire class: E / F (EN13501)
- External dimensions: W x H x D = 44 x 46 x 27 cm (*incl. door handles*)
 Internal dimensions: W x H x D = 33 x 38 x 20 cm
- Weight: Gross weight incl. packaging 18.5 kg.
-

MAINTENANCE / CLEANING:

We recommend cleaning the outer surface of the cabinet at least 1-2 times per year.
 Wipe using detergent and a damp cloth or brush and rinse with clean water. Finally wipe down the cabinet with a dry cloth.

If the surface of the cabinet becomes tarnished, this will be due to deposits which can easily be cleaned off.
 For difficult deposits, use a scouring pad. **Never use wire wool.**