


SCREENi®

Viewing solution for Broncoflex®

- High resolution screen
- Compact and mobile system
- Easy to use and intuitive
- Examination management interface
- Export photos and videos via USB



SCREENi® - Viewing solution for Broncoflex®



HIGH RESOLUTION AND MOBILITY

High resolution screen 1280x800 with a bright touch panel offering optimal colour reproduction and contrast.

Compact and mobile system easily transportable over examination sites, the display may also be exported to an already existing medical monitor.

Simple to use and intuitive, starts in seconds and immediately displays the video from the connected device. Videos and photos can be recorded during examination.

Examination Management Interface making it possible to view photos and videos of performed exams

Export recorded photos and videos via USB for inclusion in the patient file.

ACCESSORIES



Robust mobile cart designed to connect to a SCREENi® tablet and fitted with provisions on which to hang Broncoflex® bronchoscopes in their sterile packaging.

Réf. 006.0003



Universal clamp for attaching SCREENi® to rails and tubes.

Réf. 006.0020

TECHNICAL SPECIFICATIONS

Electrical	Power supply	External mains AC-DC Adapter 100 - 240V - 0,5 A 50 - 60 Hz
	Type of protection against electrocution	Class 2
	Battery type	Lithium-Ion battery
	Autonomy	3h30
	Charging time	2h
Touch screen	Screen size	10 inches (21.6 x 13.5 cm)
	Resolution	1280 x 800
Miscellaneous	IP Code	IP 30
	Dimensions (L x H x D)	310 x 220 x 770 mm
	Weight	2 Kg
	Accessories delivered with SCREENi®	Charger
	Processor	Embedded smart system
	Physical memory size	16 GB
	Input connector	Proprietary video connector 2 x USB Type A
Output connector	HDMI	
Compliance	MDD classification	Class IIa
Ordering information	SCREENi® 10" (including power cord and battery charger)	005.0008