

**S-CUT<sup>®</sup>**

# REVOLUTIONARY CUTTING TOOL



**A lifesaver in accidents and traumas**



**s-cut.se • esequipment.se**



For those who are working to save lives:

# A tool you can trust every time!

*The simple and ingenious design of the  
S-CUT allows you to trust it every time.*



*Easy to clean and maintain.  
The cutting blade can be  
removed without use of any  
tool and your S-CUT is also  
very easy to disassemble for  
cleaning.*

*Saves precious time. S-CUT helps you to immediately get access to do life saving efforts. Precious time is saved in every trauma situation.*

#### **Cuts through**

Combat vests, motorcycle suits, diving suits, several layers of fabric, wet or/and oily materials and more.

#### **You don't have to be strong**

The S-CUT tool is designed to be used by all kind of medical rescue staff or first responder, in spite of individual strength. You can cut towards or away from yourself.

#### **Safe for you and your patient**

S-CUT is designed to protect yourself as well as the patient from injury. Can't be used as a weapon if someone takes it from you.

#### **Works in all weather and temperatures**

Your S-CUT works as good in a hot desert as in an arctic environment. It is neither sensitive to heat, cold or humidity.

#### **Requires little to no education**

Since it is easy to operate, you can start using your S-CUT immediately after reading the instructions.

#### **With or without gloves**

S-CUT is easy to use even if you are wearing heavy gloves.

#### **Left- or righthanded**

Designed for ergonomics and optimum grip. As easy to use even if you are right- or lefthanded.





# S-CUT is a life saver

*S-CUT is an optimal tool when it comes to quickly removing clothing to start lifesaving measures. Compared with scissors and knives, the S-CUTs considerably shorten the time for cutting clothes in trauma situations, cardiac arrests, chemical accidents and hypothermia.*



Oil recovery at Tjörn 2011\*



*S-CUT is safe to use. Despite the extreme cutting ability there is no risk of cut injuries, neither of patient or rescue staff.*



*The 501, 601 and 701 models are intended for medical use and have high standards. Ambulance- and rescue teams all over the world are using them.*



*Originally designed for emergency-, ambulance- and rescue teams. S-CUT has now found a place within police- and, military, correctional- and decontamination personnel.*

\*At the major oil disaster at the Swedish island Tjörn in 2011, S-CUT was used inter alia to remove clothes from the rescue teams after every shift. The clothes were soaked in thick, sticky oil and by cutting them up before removing the staff spared unnecessary skin contact with the oil.




In traffic accidents S-CUT is used to rapidly cut safety belts and inflatable curtains hanging in the way. Clothing is cut to properly put the neck collar in place.

Photo: West Coast Rescue Team, 2:nd best team in Rescue World Championship, always use S-CUT in action.



*”A single minute gained can contribute significantly of a patient’s chance to survival.”*

John Herlitz, Professor of Cardiology at Sahlgrenska Universitetssjukhuset in Gothenburg



Despite it's about cutting clothing to start lifesaving measures, or getting a broken leg out of a shoe or boot without causing the patient unnecessary pain, S-CUT is a useful and competent choice.



...POL-  
...f-  
...has  
...the  
...S-CUT.  
...their  
...S-CUT  
...ster,  
...hable,  
...side of  
...armor.



# Shortens the time to action

■ Studies show that the efficiency of S-CUT spare several minutes in emergency situations

**THE POLICE FORCE** use S-CUT to help those who are in need of rapid assistance, and to gain access to a patient's chest in the case of cardiac arrest. Moreover, if a policeman is injured and needs immediate treatment, S-CUT is capable of cutting easily through the heavy outer garments of a police uniform.

**CORRECTIONAL.** Detainees and prisoners that hurt themselves are in need of quick first aid. It can involve penetrating injuries or hanging attempt. S-CUT is far more efficient compared to scissors and knives and it can't be used as a weapon to cut or stab the staff.

**RESCUE TEAMS** have great results with S-CUT's since they can be used wearing heavy gloves. In extrications and in animal rescues, straps, ropes and belts can be easily cut off. In traffic accidents it is fast and easy to cut off seat belts and inflatable curtains. The 601 model is especially useful in e.g. car crashes, where it can cut open the fabric of the back of a front seat and clear the way for the use of hydraulic scissors.

**WITHIN THE MILITARY DEFENCE** S-CUT's ergonomic design is suitable under tactically difficult conditions when it is necessary to keep a low position. In those circumstances it is an ad-

vantage that S-CUT can be used with left as well as right hand. The 601 model has an additional blade with slitting function that is activated by a slight pressure. That function eliminates the need for a free edge to start cutting and makes it possible to quickly treat the injured for e.g. gunshot wounds without exposing them to cold.

**AMBULANCES AND EMERGENCY ROOMS.** All trauma patients, cardiac arrest patients and critically ill patients require clothes to safely and quickly be removed for efficient triage and treatment. Patients with hypothermia, sometimes with wet clothes, must be handled with great care but still effectively. S-CUT's makes it possible to remove clothing in a few seconds without risk to the patient. In the ambulance, when a fast iv-cannula and/or vital parameters are required, the whole sleeve is cut from the shoulder. The patient is not exposed to cold and continued treatment can be performed in the heated ambulance.

**DECONTAMINATION.** Before decontamination can be done clothes must be removed, on laying as well as on standing people. With an S-CUT the clothes are cut off to avoid pulling them over the head. The tool works on wet clothes as well as on clothes smeared with thick oil.

# Voices about S-CUT

■ Our best reviews come from professionals who use S-CUT in varied situations

*Our products constantly find new groups of users. Departments within emergency care, military and police already know the benefit of having an S-CUT available. We are proud to refer to the people who require high standards of functioning equipment. Each and every one of them is well aware that a patient's life and health after accidents and injuries depends on how quick clothes can be removed and life saving efforts can be performed.*

”

S-CUT is an important part of the equipment I use to take care of injured soldiers. I have used it for years and find it totally reliable. S-CUT's ability to cut through the soldiers equipment, including the body armor, is unique.

Mats, Reconnaissance soldier/nurse



”

The best thing about S-CUT is that it is so quick and easy to use. You can use it in all types of situations, from the most routine to the most extreme.

Martin Sjöholm, health care manager,  
Swedish Rescue services international relief efforts



”

Of all tools we have tested in accidents and traumas, S-CUT have time after time proved that it is the best. For example, we save up to five minutes to get the leather suit of a motorcyclist.

Sari Sirilää, Nurse/ATLS coordinator,  
Sahlgrenska University Hospital Göteborg.



**More reviews on [esequipment.se](http://esequipment.se)**



# More products from ES Equipment



## S-CUT 501

Our largest model, designed for really tough conditions. Medical device with high demands on precision and performance. Length 24 cm, weight 260 g. Anodized machine milled aluminum.



## S-CUT QE

An allround cutting device for those who need to cut through different materials, almost everything but metal. With amazing cutting capability you can easily cut heavy leather belt, tarpaulin or ground cloth. Length 18 cm, weight 106 g.



## S-CUT 701

Our smallest model that also is designed for really tough conditions. Medical device with high demands on precision and performance. Length 18 cm, weight 170 g. Anodized machine milled aluminum.



## Holster

Made of durable fabric that can handle the toughest conditions. Washable. Belt clip on the back. For all S-CUT models.



## S-CUT 601

As 701 but completed with a extendable blade for making hole to start the cutting. Medical device with high demands on precision and performance. Length 18 cm, weight 170 g. Anodized machine milled aluminum.



## Wall mount

Keep your S-CUT readily available in this practical wall mount! You know where you have it, and in one single grab your S-CUT is ready for emergency situations. For all S-CUT models.



## S-CUT XC-E – available in four colors

Adopted for situations where horses get stuck and panic. A fast and safe way to relieve the horse without injure it or yourself. Has a bigger gap to fit thick halter straps and saddle girths etc. Length 21 cm, weight 105 g. Powder coated water-cut aluminium.



## Cutting blade and plastic support

Blister pack with 1 cutting blade and 1 plastic support.

**For S-CUT XC-E.**



## ACP

4-pack of Blades and Hygiene Supports.

**For model 501, 601 and 701.**