

Pediatric Solutions

ZOLL®



Optimal Support for
Pediatric Resuscitation

Because Children Are Not Just Little Adults

Too often, clinicians who treat children need to compromise and adapt devices designed for adults because there just isn't a pediatric version available. That is not the case with the ZOLL® R Series®, which has unique pediatric capability built into the defibrillator to optimize resuscitation of children. With smart, simple, and safe pediatric support, whether you need a solution for a single pediatric department or an entire pediatric hospital, the ZOLL R Series offers the most extensive pediatric capabilities available in a professional defibrillator.

The pediatric electrodes can be used for monitoring, pacing, defibrillation, cardioversion, and CPR.



R Series Supports Pediatric Resuscitation

At ZOLL, we recognize that the resuscitation needs of children are different from adults and have provided a range of capabilities to support the special needs of the youngest patients.



CPR Support Optimizes Perfusion and Reduces Interruptions

With the vast majority of pediatric arrests due to respiratory failure, the importance of high-quality CPR cannot be minimized. That's why the R Series is the first device to provide a pediatric electrode with a built-in CPR sensor. It is difficult to judge how hard one is pushing when performing CPR on a child, so the CPR Dashboard™ reports the actual depth and rate of compressions delivered. A CPR timer gauges the CPR periods for optimal ventilation, a metronome activates during manual mode when the rate falls below 100 compressions per minute, and an idle timer keeps staff aware of interruptions that compromise perfusion. All of the CPR information is contained in a single location at the top of the screen for rapid assessment.

See-Thru CPR® is a proprietary filter that allows clinicians to see if an organized rhythm is developing during CPR, helping to reduce the duration of interruptions to perfusion and speeding time to shock when indicated.



Pediatric AED Algorithm Adjusts for Pediatric Rhythms

In the event of an arrhythmic arrest, as soon as the pediatric electrodes are placed, the R Series automatically adjusts its AED analysis algorithm to pediatric parameters in order to reduce the likelihood that a compensating rhythm will be terminated, enhancing the likelihood that unique pediatric ECG morphology will be identified when a shock or cardioversion is indicated.

Constant Current Pacing for Bradycardia

ZOLL's constant current pacer has been proven to deliver more capture with less current than competing devices, ensuring support for the bradycardia often seen in children in distress¹.

Optional Monitoring Devices Have Pediatric Sensors

Neonatal and pediatric cuffs for non-invasive blood pressure measurement and infant and pediatric SpO₂ sensors provide important monitoring information for evaluating potential patient deterioration, they can be equally useful during transport or resuscitation. Optional mainstream EtCO₂ monitoring using the Respirationics Capnostat 5 Sensor with pediatric airway adapters help confirm and monitor proper endotracheal tube placement and alerts to any potential dislodgement. It also provides valuable information on the progress of resuscitation.



Focus on Simplicity Minimizes Errors and Confusion

The R Series keeps things simple for clinicians. The pediatric paddles are stored inside the adult paddles. And the R Series uses the same OneStep® cable to connect both adult and pediatric electrodes as well as internal spoons.

Safe Operation Is a Priority



When the pediatric electrodes are attached to the R Series, the shock energy is automatically lowered to a **starting**

dose of 50 joules. For units

like the NICU, where even 50

joules may be too high, the R series can be configured to automatically decrement to whatever dose is desired, **even as low as 1 joule.** By automatically decreasing the starting dose, the R Series, minimizes the chance of incorrect dosing. This is different from other devices that attenuate the energy by adding resistance to the circuit, which impacts impedance measurement and reports inaccurate dosage.

Automatic Testing Assures Defibrillator Readiness

Pediatric resuscitation is infrequent, but devices must always be ready. The R Series automatically performs a daily 100-point self test, with the electrodes attached. The test will even alert to an expired electrode. If the unit is compromised, the cause is displayed on the screen. If desired, the results of the self-test can be wirelessly transmitted to clinical engineering to assist in troubleshooting.

Don't Compromise When It Comes to Pediatrics...Rely on ZOLL

¹ZOLL PM, et al. *Circulation*. May;71(5):937-44.

ZOLL Medical Corporation, an Asahi Kasei Group company, develops and markets medical devices and software solutions that help advance emergency care and save lives, while increasing clinical and operational efficiencies. With products for defibrillation and monitoring, circulation and CPR feedback, data management, fluid resuscitation, and therapeutic temperature management, ZOLL provides a comprehensive set of technologies that help clinicians, EMS and fire professionals, and lay rescuers treat victims needing resuscitation and critical care. For more information, visit www.zoll.com.

The Asahi Kasei Group is a diversified group of companies led by holding company Asahi Kasei Corp., with operations in the chemicals and fibers, homes and construction materials, electronics, and health care business sectors. Its health care operations include devices and systems for critical care, dialysis, therapeutic apheresis, transfusion, and manufacture of biotherapeutics, as well as pharmaceuticals, diagnostic reagents, and nutritional products. With more than 25,000 employees around the world, the Asahi Kasei Group serves customers in more than 100 countries. For more information, visit www.asahi-kasei.co.jp/asahi/en/.

**ADVANCING
RESUSCITATION.
TODAY.®**

ZOLL Medical Corporation
Worldwide Headquarters
269 Mill Road
Chelmsford, MA 01824
978-421-9655
800-348-9011
www.zoll.com

For subsidiary addresses and fax numbers, as well as other global locations, please go to www.zoll.com/contacts.