


# Nuflow®

## nasal cannula



Available sizes:

	1 Nuflow® nasal cannula (small)	2 Nuflow® nasal cannula (medium)	3 Nuflow® nasal cannula (large)	4 Nuflow® nasal cannula (x-large)
REF	1320	1330	1340	1350
low limits <sup>1</sup>	10,0 l/min	14 l/min	23 l/min	27 l/min
Outer diameter prong	2,3 mm	2,7 mm	3,2 mm	3,7 mm
Septum distance	3,7 mm	4,3 mm	4,8 mm	5,3 mm
Patient group <sup>2</sup>	Premature	Premature, Newborns	Infants	Children
Packing Unit	10 pieces	10 pieces	10 pieces	10 pieces
Fixation patches (Spare part)	1380 small	1380 small	1381 large	1381 large

### Connectors

Tube connection: OD 10 mm or OD 15 mm available

1 - See user manual for details

2 - Depends on size of patient nostril



## Material

Cannula Silicone, Zylar, TPE, PP,  
Fixation patch Hydrocolloid  
- free of PVC, DEHP, natural rubber latex

---

## Fixation patches (spare part)

REF 1380 Pair of Fixation-patches (small), (spare part for REF 1320/ REF 1330)  
REF 1381 Pair of Fixation-patches (large), (spare part for REF 1340/ REF 1350)

---

## Starter Sets

REF 1393 Nuflow® starter set s/m (5x Nuflow® nasal cannula (small),  
5x Nuflow® nasal cannula (medium))  
REF 1394 Nuflow® starter set l/xl (5x Nuflow® nasal cannula (large),  
5x Nuflow® nasal cannula (x-large))

---

## Combinable devices

medinSINDI®, medinCNO®mini, medin-NC3®, medinBLENDER

---

## Regulations

Classification	EU IIa
Country of origin	Germany
Intended use	single patient use, non invasive, max. 7 days
Manufacturing surrounding	clean room
Disposal	The Nuflow® nasal cannula and its packaging can be discarded with regular household waste.



Manufacturer:

medin Medical Innovations GmbH

Adam-Geisler-Straße 1

82140 Olching, Germany

+49 (0)8142 448 46 -0

info@medin-medical.com

www.medin-medical.com



CE0483

REV03  
17.06.2020

Specifications are subject to change without notice. Some features are options. Not all features are available in all markets.  
For all proprietary trademarks (®) and third-party trademarks (®) used by medin Medical Innovations GmbH see  
www.medin-medical.com/trademarks. © 2020 medin Medical Innovations GmbH. All rights reserved.

Nuflow