



High quality products for rescue

Our **NEW** premium product **S-CUT 501-RTL**, an improved version of the old 501, is an optimal tool when it comes to start rescue measures. Compared with scissors and knives, the **S-CUT** considerably shorten the time for cutting clothes in rescue situations, cardiac arrests, chemical accidents and hypothermia.

S-CUT is safe to use. Despite the extreme cutting ability the risk of cut injuries is almost non existent, neither of patient nor rescue staff.

The **501-RTL** model is intended for rescue teams and therefor produced with premium quality. Ambulance-, rescue teams

and military units all over the world are using it and more and more organisations discover that **S-GUT** is unbeatable.

S-CUT has also found a place within police, correctional- and decontamination personnel.

A tool you can trust every time

The simple and ingenious design allows you to trust it every time

Voices about S-CUT

- The best thing about **S-CUT** is that it is so quick and easy to use. You can use it in all types of situations, from the most routine to the most extreme.

 Martin Sjöholm, health care manager, Swedish Rescue services international relief efforts
- Of all tools we have tested in accidents and traumas, **S-CUT** have time after time proved that it is the best. For example, we save up to five minutes to get the leather suit of a motorcyclist.

 Sari Sirilää, Nurse/ATLS coordinator, Sahlgrenska University Hospital Göteborg
- **S-CUT** is an important part of the equipment I use to take care of injured soldiers. I have used it for years and find it totally reliable. **S-CUT's** ability to cut through the sodiers equipment, including the body armor, is uniqe.

 Mats, Operator

S-CUT 501-RTL

Weight: 0,260 kg. Length: 240 mm.

Material: Anodized machine milled

aluminum.

Accessories: Wall mount and holster.



S-CUT is developed, produced and marketed by ES Equipment AB







More from us

701 – Rescue device for high demands

The S-CUT 701 model is popular by ambulance personnel due to the smaller size. Like 501-RTL it is a rescue device for high demands on precision and performance and have the same exeptional cutting abilities.



S-CUT 701

Weight: 0,170 kg. Length: 180 mm.

Material: Anodized machine milled

aluminum.

Accessories: Wall mount and holster.

ACP-RTL

A 4-pack of Blades and Hygiene Supports for 501-RTL.

Make sure you always have these parts available!



ACP

A 4-pack of Blades and Hygiene Supports for 501, 601 and 701.

Make sure you always have these parts available!



