




# STATPACKS

## THE G3 BREATHER IS AN ALL-INCLUSIVE EMS AIRWAY MGMT PACK FOR AIRWAY CONTROL AND O2 THERAPY.

The Breather is a versatile EMS airway management pack that can be used in many different departments and settings.



-  **ACLS** (Advanced Cardiovascular Life Support)
-  **AED** (Automated Extrenal Defibrillator)
-  **ALS** (Advanced Life Support)
-  **BVM** (Bag Valve Mask)
-  **PAL** (Pediatric Advanced Life Support)
-  **BLS** (Basic Life Support)
-  **O2** (Accommodates D or Jumbo D Oxygen Cylinder)



GO ANYWHERE. DO ANYTHING.



## THE G3 BREATHER MEETS LOCAL PROTOCOLS FOR OXYGEN AND AIRWAY REQUIREMENTS.

StatPacks is a top-notch EMS bag system that is built with the men and women of the Fire/EMS community in mind. The G3 Responder and the G3 BackUp are on every front-line apparatus in my department, and will be for many years to come!

**Michael Hadden**  
Asst. Division Chief EMS/Training - Tallahassee Fire Department



RED



GREEN



BLUE



BLACK

[OVERVIEW](#) [SPECIFICATIONS](#) [BUY](#)

### G3 Breather

- Contents NOT included
- Unique design allows access to the main compartment from the top or front of this EMS airway management Backpack
- Specifically designed to hold StatPacks Generation 3 Emergency Medical Services Kits; including all StatPacks® First Aid QuickRoll Intubation, Remedy, and Circulatory kits
- QuickZip access to front compartment for quick access to EMS airway management Emergency Medical Equipment
- Tarpaulin bottom panel for protection from abrasion and the elements
- Clear, durable urethane windows, mesh pockets, and elastic holsters keep contents in place and well-organized
- Added more durable reflective material that has a protective outer mesh cover which shields it from wear and tear. Helps with visibility in nighttime situations (Not included in Tactical Black)
- Accommodates removable O2 Module (sold separately)
- Foam molded shoulder straps that allow for comfortable transport of most emergency medical equipment
- Internal Velcro® lining strips to secure modules or other Velcro® compatible items to secure ambulance equipment
- LED light Attachment point (lights sold separately)
- Unique foam lined design keeps all Emergency Medical Equipment inside protected from breakage or the elements
- Side pockets for easy access to vital sign tools, bandages, and smaller ambulance equipment
- A quick stash pocket for trash and extra supplies located in the outside front
- ID window to label any way you want to identify between other StatPacks EMS Backpack's
- Main fabric is a non-toxic PVC-free tarpaulin. An environmentally friendly fabric featuring easy-clean, pathogen and water resistant properties

### Size and Weight

- Height:** 21 inches
- Width:** 19 inches
- Depth:** 7 inches
- Weight:** 4.5 lbs

### Available Colors

Red, Green, Blue, Tactical Black

### SKU

- SKU: G35008RE [red]
- SKU: G35008GN [green]
- SKU: G35008BU [blue]
- SKU: G35008TK [tactical black]