

ZOLL AED defibrillator accessories

CONTENT

- 2 | Electrodes
- 3 | Batteries
- 4 | Storage
- 5 | AED Cabinets
- 6 | Training






8000.0839



8900.0003



8900.0002.70



8900.0400 / 8900.0402



8900.0800.01



8900.0801.01 / 8900.0802.01



8900.0810.01



8900.0807.01

Electrodes

			PLUS	PRO	AED 3
8000.0839	AED Pro ECG cable. European. 3-lead ECG cable	1 pc		X	
8900.0003	ECG electrodes for AED Pro	Box of 30 pcs		X	
8900.0002.70	CPR Uni-Padz II Electrode. Universal with feedback. Pediatric/adult. 5 years shelf life	1 pc			X
8900.0400	CPR Stat-Padz Electrode. Multi-Function. Adult. 2 years shelf life	Box of 8 pcs	X	X	X
8900.0402	CPR Stat-Padz Electrode. Multi-Function. Adult. 2 years shelf life	1 pc	X	X	X
8900.0800.01	CPR-D Padz One-Piece Electrode Pad with feedback. Adult. Incl. CPR-D Accessory kit. 5 years shelf life	1 pc	X	X	X
8900.0801.01	Stat-Padz II Electrode. Multi-Function. Adult. 2 years shelf life	1 pc	X	X	X
8900.0802.01	Stat-Padz II Electrode. Multi-Function. Adult. 2 years shelf life	Box of 12 pcs	X	X	X
8900.0810.01	Pedi-Padz II Electrode. Multi-Function. Pediatric. 1 < 8 year or < 25 kg. 2 years shelf life	1 pc	X	X	X
8900.0807.01	CPR-D Accessory Kit. CPR Barrier Mask, scissor, gloves, prep razor, towel and a moist towelette	1 pc	X	X	





8000.0006.96



8000.0807.01



8000.0860.01



8019.0535.01



8200.000100.XX

Batteries

			PLUS	PRO	AED 3
8000.0006.96	AED 3 battery	1 pc			X
8000.0807.01	AED Plus battery. Type 123A lithium-mangan batteries. 5 years shelf life in standby	Box of 10 stk	X		
8000.0860.01	AED Pro lithium battery. Non-rechargeable. 5 years shelf life	1 pc		X	
8019.0535.01	AED Pro lithium battery. SurePower. Rechargeable. 5 years shelf life	1 pc		X	
8200.0001.00XX	SurePower Single Bay Charger for 8019.0535.01	1 pc		X	



8000.0012.50



8000.0012.53



8000.0012.54



8000.0012.55



8000.0012.66



8000.0802.01



8000.0809



8000.0810.01



8000.0832.01



8000.0875.32

Storage

			PLUS	PRO	AED 3
8000.0012.50	AED 3 Carry case	1 pc			X
8000.0012.53	AED 3 Pelican case. Small	1 pc			X
8000.0012.54	AED 3 Pelican case. Large	1 pc			X
8000.0012.55	AED 3 Mounting bracket	1 pc			X
8000.0012.66	AED 3 Mounting bracket. Space for a AED 3 with carry case	1 pc			X
8000.0802.01	AED Plus carry case	1 pc	X		
8000.0809	AED Plus Mounting bracket	1 pc	X		
8000.0810.01	AED Pro Soft carry case. Nylon. Standard	1 pc		X	
8000.0832.01	AED Pro Carry case. Vinyl. Spare battery compartment (battery, electrode and 3-lead ECG cable)	1 pc		X	
8000.0875.32	AED Pro Pelican case	1 pc		X	





CX-OD-GN-XX



CX-OD-WE-XX



CX-OD-GY-XX



CX-ID-B-XX



CA HSS100



CA HSS100P



CA HSS100GSM-EXP



CA HSS101GSM-EXP



8000.0855



8000.0012.56



Aivia100-XX



Aivia200-XX

AED Cabinets

			PLUS	PRO	AED 3
CX-OD-GN-XX	Cabix Outdoor AED Cabinet. Outdoor. Monitoring. Alarm. Light. Heat. Climate friendly. Power saving. Green	1 pc	X	X	X
CX-OD-WE-XX	Cabix Outdoor AED Cabinet. Outdoor. Monitoring. Alarm. Light. Heat. Climate friendly. Power saving. White	1 pc	X	X	X
CX-OD-GY-XX	Cabix Outdoor AED Cabinet. Outdoor. Monitoring. Alarm. Light. Heat. Climate friendly. Power saving. Grey	1 pc	X	X	X
CX-ID-B-XX	Cabix Indoor AED Cabinet. Indoor. Alarm. Battery indicator. Power saving	1 pc	X	X	X
CA HSS100	HSS100 AED Cabinet. Outdoor. Alarm. Light. Heat. Ventilation	1 pc	X	X	X
CA HSS100P	HSS100 AED Cabinet. Outdoor. Alarm. Light. Heat. Ventilation. Extreme cold -40 °C	1 pc	X	X	X
CA HSS100GSM-EXP	HSS100GSM AED Cabinet. Outdoor. Alarm. Light. Heat. Ventilation. GSM/SMS communication	1 pc	X	X	X
CA HSS101GSM-EXP	HSS101GSM AED Cabinet. Outdoor. Alarm. Light. Heat. Ventilation. GSM/SMS communication. AED sensor	1 pc	X	X	X
CA HSSST2	Stand for CA HSS100 AED Cabinet. Stainless steel	1 pc	X	X	X
CA HSS STRS	Stand for CA HSS101GSM AED Cabinet. Stainless steel	1 pc	X	X	X
8000.0855	ZOLL AED Plus Wall cabinet. Indoor. Alarm. White metal wall cabinet (W: 44,5 x H: 44,5 x D: 23 cm)	1 pc	X		
8000.0012.56	ZOLL AED 3 Wall cabinet. Indoor. Alarm. White metal wall cabinet (W: 44,5 x H: 44,5 x D: 18 cm)	1 pc			X
Aivia100-XX	Aivia 100 AED Cabinet. Indoor. Alarm. Battery.	1 pc	X	X	X
Aivia200-XX	Aivia 200 AED Cabinet. Outdoor. Alarm. Light. Heat	1 pc	X	X	X



185.10050



8000.0009.25



8900.0002.74



8900.0002.75



8028.0000.01XX



8028.0000.10



8028.0000.12



8028.0000.13



8008.0050.XX



8900.0190



8900.0195



8900.5007



8900.0803.01



8900.0804.01

Training

			PLUS	Trainer 2	AED 3	AED 3 Trainer	PRO
185.10050	Shocklink System (requires adapter 185-50450)	1 pc	X		X		X
8000.0009.25	AED Simulator, English	1 pc	X		X		X
8900.0002.74	CPR Uni-Padz II Training electrode. Y connector for simulator connection	1 pc			X		
8900.0002.75	CPR Uni-Padz II Training electrode. Replacement pads for 8900.0002.74	Box of 8 pcs			X		
8028.0000.01XX	ZOLL AED 3 Trainer unit	1 pc				X	
8028.0000.10	CPR Uni-Padz Training electrode harness	1 pc				X	
8028.0000.12	CPR Uni-Padz Training electrode replacement adhesive gels	Box of 5 pcs				X	
8028.0000.13	CPR Uni-padz Training electrode replacement liners for 8028.0000.12	1 pc				X	
8008.0050.XX	ZOLL AED Plus Trainer 2 unit	1 pc		X			
8900.0190	CPR Stat-Padz Training electrode. Y connector for simulator connection	1 pc	X				X
8900.0195	CPR Stat-Padz Training electrode. Replacement pads for 8900.0190	Box of 8 pcs	X				X
8900.5007	CPR-D Padz Training electrode. Y connector for simulator connection. Gel pads	1 pc	X				X
8900.0803.01	Replacement Gel Pads for CPR-D Padz training electrode (8900.0804.01 and 8900.5007)	Box of 5 pcs	X	X			X
8900.0804.01	CPR-D Padz Training electrode	1 pc		X			
8900.0805.01	Stat-Padz II Training electrode. Adult. 2 years shelf life	Box of 6 pcs		X			
8900.0008.6101	Pedi-Padz II Training electrode. Pediatric. 2 years shelf life	Box of 6 pcs		X			

