


# HELPING YOU PREVENT LOWER RESPIRATORY TRACT INFECTIONS



# SOLUTIONS FOR RESPIRATORY HEALTH

Avanos is committed to help you reduce lower respiratory tract infections by providing effective clinical solutions to help with the removal and prevention of aspiration and colonisation in the oropharynx.

Avanos's Respiratory Health solutions bundle includes the BALLARD<sup>®</sup> Closed Suction Systems, MICROCUFF<sup>®</sup> Endotracheal Tubes, BALLARD<sup>®</sup> Oral Care Systems and BALLARD<sup>®</sup> Mini-BAL Sampling Catheter.



# TABLE OF CONTENTS

<b>CLOSED SUCTION SYSTEMS - ADULTS</b>	5
<b>BALLARD* TURBO-CLEANING CLOSED SUCTION SYSTEMS</b>	5
BALLARD* Turbo-Cleaning Closed Suction System - Double Swivel Elbows	5
BALLARD* Turbo-Cleaning Closed Suction System - WET PAK* Double Swivel Elbows	7
BALLARD* Turbo-Cleaning Closed Suction System - T-Pieces	8
BALLARD* Turbo-Cleaning Closed Suction System - WET PAK* T-Pieces	9
<b>BALLARD* MULTI-ACCESS PORT CLOSED SUCTION SYSTEMS</b>	10
BALLARD* Multi-Access Port Closed Suction System - Double Swivel Elbows	10
BALLARD* Multi-Access Port Closed Suction System - WET PAK* Double Swivel Elbows	12
BALLARD* Multi-Access Port Closed Suction Replacement Catheter	14
<b>BALLARD* CLOSED SUCTION SYSTEMS</b>	15
BALLARD* Closed Suction System - T-Pieces	15
BALLARD* Closed Suction System - Double Swivel Elbows	16
BALLARD* Closed Suction System - WET PAK* Double Swivel Elbows	17
<b>BALLARD* CLOSED SUCTION SYSTEMS - ACCESSORIES</b>	18
BALLARD* Multi-Access Port Closed Suction System - Accessories	18
BALLARD* Closed Suction System - Accessories	18
BALLARD* Closed Ventilator Circuit Drain	19
<b>ENDOTRACHEAL TUBES - ADULTS</b>	20
MICROCUFF* Subglottic Suctioning Endotracheal Tubes	20
MICROCUFF* Endotracheal Tubes	20
<b>ORAL CARE SOLUTIONS - ADULTS</b>	21
<b>BALLARD* ORAL CARE KITS</b>	21
BALLARD* Oral Care Advanced Kit	21
BALLARD* Oral Care Moderate Kit	21
BALLARD* 24-Hour Oral Care Kit (q4)	22
BALLARD* Mini Oral Care Kit	23
<b>BALLARD* ORAL CARE SUCTION KITS - COMPONENTS &amp; SUCTION KITS</b>	24
BALLARD* Oral Care Suction Kits	24
BALLARD* Oral Care Suction Components	25
BALLARD* Oral Care Swabs	26
BALLARD* Mini Oral Care Swabs	26
<b>DIAGNOSTIC CATHETERS - ADULTS</b>	27
<b>BALLARD* MINI-BAL SAMPLING CATHETER</b>	27

<b>CLOSED SUCTION SYSTEMS - NEONATES/PAEDIATRICS</b>	28
BALLARD* Closed Suction System for Neonates/Paediatrics - Elbows	28
BALLARD* Closed Suction System for Neonates/Paediatrics - "Y"	29
BALLARD* Multi-Access Catheter for Neonates/Paediatrics	30
BALLARD* Closed Suction System for Neonates/Paediatrics - Accessories	30
<b>ENDOTRACHEAL TUBES - PAEDIATRICS</b>	31
MICROCUFF* Endotracheal Tube for Paediatric, Oral/Nasal Magill	31
MICROCUFF* Endotracheal Tube for Paediatric, Oral Curved	32
<b>ORAL CARE SOLUTIONS - NEONATES/PAEDIATRICS</b>	33
BALLARD* Oral Care Suction Catheters for Neonates/Paediatrics	33
<b>INDEX BY CODE</b>	34
<b>NOTES</b>	35


# CLOSED SUCTION SYSTEMS - ADULTS

## BALLARD\* TURBO-CLEANING CLOSED SUCTION SYSTEMS

**Indication:** Aspiration catheter for removal of tracheal secretions of intubated / tracheotomized adult patients respecting the closed ventilation circuit.

### BALLARD\* TURBO-CLEANING CLOSED SUCTION SYSTEM - DOUBLE SWIVEL ELBOWS

- BALLARD\* Technology
- Built-in Cleaning Chamber
- Gamma Sterilised: 5 Years Shelf Life
- Available for Endotracheal and Tracheostomy Tubes
- Available with Metered Dose Inhaler (MDI) Port
- Available with dual lumen catheter
- Validated for a maximum Usage Time of 72 Hours
- Does not contain Natural Rubber Latex
- Does not contain DEHP
- Connector 15mm Female to 22mm Female



DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal MDI Port	<b>FSQ1074</b>	<b>2271018-5</b>	<b>10 Fr/3.3mm</b>	<b>54cm</b>	<b>black</b>	<b>20/case</b>
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal	<b>FSQ1076</b>	<b>227101-5</b>	<b>10 Fr./3.3mm</b>	<b>54cm</b>	<b>black</b>	<b>20/case</b>
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal	<b>FSQ475</b>	<b>22716-5</b>	<b>12 Fr/4mm</b>	<b>54cm</b>	<b>white</b>	<b>20/case</b>
<b>DOUBLE SWIVEL ELBOW</b> Tracheostomy	<b>FSQ477</b>	<b>2271603-5</b>	<b>12 Fr/4mm</b>	<b>30.5cm</b>	<b>white</b>	<b>20/case</b>
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal MDI Port	<b>FSQ1078</b>	<b>2271218-5</b>	<b>12 Fr/4mm</b>	<b>54cm</b>	<b>white</b>	<b>20/case</b>
<b>DOUBLE SWIVEL ELBOW</b> Tracheostomy MDI Port	-	<b>22712183</b>	<b>12Fr/4mm</b>	<b>30.5cm</b>	<b>white</b>	<b>20/case</b>
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal	<b>FSQ476</b>	<b>227-5</b>	<b>14 Fr/4.6mm</b>	<b>54cm</b>	<b>green</b>	<b>20/case</b>


All codes in **bold** indicate a sterile product

DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
<b>DOUBLE SWIVEL ELBOW</b>	<b>FSQ478</b>	<b>22703-5</b>	<b>14 Fr/4.6mm</b>	<b>30.5cm</b>	<b>green</b>	<b>20/case</b>
Tracheostomy						
<b>DOUBLE SWIVEL ELBOW</b>	<b>FSQ1080</b>	<b>2271418-5</b>	<b>14 Fr/4.6mm</b>	<b>54cm</b>	<b>green</b>	<b>20/case</b>
Endotracheal MDI Port						
<b>DOUBLE SWIVEL ELBOW</b>	<b>FSQ1081</b>	<b>22714183-5</b>	<b>14 Fr/4.6mm</b>	<b>30.5cm</b>	<b>green</b>	<b>20/case</b>
Tracheostomy MDI Port						
<b>DOUBLE SWIVEL ELBOW</b>	<b>FSQ1082</b>	<b>227161-5</b>	<b>16 Fr/5.3mm</b>	<b>54cm</b>	<b>orange</b>	<b>20/case</b>
Endotracheal						
<b>DOUBLE SWIVEL ELBOW</b>	<b>FSQ1083</b>	<b>2271613-5</b>	<b>16 Fr/5.3mm</b>	<b>30.5cm</b>	<b>orange</b>	<b>20/case</b>
Tracheostomy						
<b>DOUBLE SWIVEL ELBOW</b>	<b>FSQ1084</b>	<b>2271618-5</b>	<b>16 Fr/5.3mm</b>	<b>54cm</b>	<b>orange</b>	<b>20/case</b>
Endotracheal MDI Port						
<b>DOUBLE SWIVEL ELBOW</b>	<b>FSQ1085</b>	<b>22716183-5</b>	<b>16 Fr/5.3mm</b>	<b>30.5cm</b>	<b>orange</b>	<b>20/case</b>
Tracheostomy MDI Port						

All codes in **bold** indicate a sterile product

## BALLARD\* TURBO-CLEANING CLOSED SUCTION SYSTEM - WET PAK\* DOUBLE SWIVEL ELBOWS

- Contains in addition 12 Saline Vials
- Gamma Sterilised: 3 Years Shelf Life - Saline Vials: aseptically filled
- Does not contain Natural Rubber Latex nor DEHP
- Validated for a maximum Usage Time of 72 hours




DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal	<b>FSQ465</b>	<b>227166-5</b>	<b>12 Fr/4mm</b>	<b>54cm</b>	<b>white</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Tracheostomy	<b>FSQ467</b>	<b>22716036-5</b>	<b>12 Fr/4mm</b>	<b>30.5cm</b>	<b>white</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal MDI Port	<b>FSQ1000</b>	<b>22712186-5</b>	<b>12 Fr/4mm</b>	<b>54cm</b>	<b>white</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Tracheostomy MDI Port	<b>FSQ998</b>	<b>227121863-5</b>	<b>12 Fr/4mm</b>	<b>30.5cm</b>	<b>white</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal	<b>FSQ466</b>	<b>2276-5</b>	<b>14 Fr/4.6mm</b>	<b>54cm</b>	<b>green</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Tracheostomy	<b>FSQ470</b>	<b>227036-5</b>	<b>14 Fr/4.6mm</b>	<b>30.5cm</b>	<b>green</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal MDI Port	<b>FSQ1002</b>	<b>22714186-5</b>	<b>14 Fr/4.6mm</b>	<b>54cm</b>	<b>green</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Tracheostomy MDI Port	<b>FSQ1001</b>	<b>227141863-5</b>	<b>14 Fr/4.6mm</b>	<b>30.5cm</b>	<b>green</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal	<b>FSQ1088</b>	<b>2271616-5</b>	<b>16 Fr/5.3mm</b>	<b>54cm</b>	<b>orange</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal MDI Port	<b>FSQ1090</b>	<b>22716186-5</b>	<b>16 Fr/5.3mm</b>	<b>54cm</b>	<b>orange</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Tracheostomy MDI Port	<b>FSQ1091</b>	<b>227161863-5</b>	<b>16 Fr/5.3mm</b>	<b>30.5cm</b>	<b>orange</b>	<b>20/case</b>

All codes in **bold** indicate a sterile product

## BALLARD\* TURBO-CLEANING CLOSED SUCTION SYSTEM - T-PIECES

- Gamma Sterilised: 5 Years Shelf Life
- Does not contain Natural Rubber Latex nor DEHP
- Validated for a maximum Usage Time of 72 hours
- 22mm x 22mm Flex Tube Connector
- 15mm x 22mm Swivel Adapter (except 227105-5<sup>◇</sup>)




DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
<b>T-PIECE</b> Endotracheal	<b>FSQ1092</b>	<b>227105-5<sup>◇</sup></b>	<b>10 Fr/3.3mm</b>	<b>54cm</b>	<b>black</b>	<b>20/case</b>
<b>T-PIECE</b> Tracheostomy	<b>FSQ1093</b>	<b>2271053-5</b>	<b>10 Fr/3.3mm</b>	<b>30.5cm</b>	<b>black</b>	<b>20/case</b>
<b>T-PIECE</b> Tracheostomy	<b>FSQ1095</b>	<b>227135-5</b>	<b>12 Fr/4mm</b>	<b>30.5cm</b>	<b>white</b>	<b>20/case</b>
<b>T-PIECE</b> Endotracheal	-	<b>22715-5</b>	<b>12Fr/4mm</b>	<b>54cm</b>	<b>white</b>	<b>20/case</b>
<b>T-PIECE</b> Tracheostomy	-	<b>2270135</b>	<b>14Fr/4.6mm</b>	<b>30.5cm</b>	<b>green</b>	<b>20/case</b>

All codes in **bold** indicate a sterile product



## BALLARD\* TURBO-CLEANING CLOSED SUCTION SYSTEM - WET PAK\* T-PIECES

- Contains in addition 12 Saline Vials
- Gamma Sterilised: 3 Years Shelf Life - Saline Vials: aseptically filled
- Validated for a maximum Usage Time of 72 hours
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
<b>WET PAK* T-PIECE</b> Endotracheal	<b>FSQ1100</b>	<b>2271516-5</b>	<b>12 Fr/4mm</b>	<b>54cm</b>	<b>white</b>	<b>20/case</b>
<b>WET PAK* T-PIECE</b> Endotracheal	<b>FSQ1102</b>	<b>2270156-5</b>	<b>14 Fr/4.6mm</b>	<b>54cm</b>	<b>green</b>	<b>20/case</b>
<b>WET PAK* T-PIECE</b> Tracheostomy	<b>FSQ1103</b>	<b>22701356-5</b>	<b>14 Fr/4.6mm</b>	<b>30.5cm</b>	<b>green</b>	<b>20/case</b>
<b>WET PAK* T-PIECE</b> Tracheostomy	-	<b>22012356-5</b>	<b>14 Fr/4.6mm</b>	<b>30.5cm</b>	<b>green</b>	<b>20/case</b>


All codes in **bold** indicate a sterile product

# BALLARD\* MULTI-ACCESS PORT CLOSED SUCTION SYSTEMS

**Indication:** Aspiration catheter for intubated or tracheotomized adult patients who are on mechanical ventilation allowing doctors to perform alternate respiratory treatments without breaching the circuit.

## BALLARD\* MULTI-ACCESS PORT CLOSED SUCTION SYSTEM - DOUBLE SWIVEL ELBOWS

- BALLARD\* Technology
- Built-in Cleaning Chamber
- Gamma Sterilised: 5 Years Shelf Life
- Rotational Manifold with Multiple Ports (Suctioning Port & Alternate Therapy for Bronchoscopy or BALLARD\* Mini-BAL sampling catheter)
- Tethered Sleeve
- Tethered Cap
- Available for Endotracheal and Tracheostomy Tubes
- Available with Metered Dose Inhaler (MDI) Port
- Suction Catheter: Validated for a maximum Usage Time of 72 Hours
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal	-	<b>24010-5</b>	10 Fr/3.3mm	58.5cm	black	20/case
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal	-	<b>240121-5</b>	12 Fr/4mm	58.5cm	white	20/case
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal MDI Port	-	<b>240128-5</b>	12 Fr/4mm	58.5cm	white	20/case
<b>DOUBLE SWIVEL ELBOW</b> Tracheostomy MDI Port	-	<b>24001238-5</b>	12 Fr/4mm	35cm	white	20/case
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal	-	<b>24014-5</b>	14 Fr/4.6mm	58.5cm	green	20/case
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal MDI Port	-	<b>240148-5</b>	14 Fr/4.6mm	58.5cm	green	20/case


All codes in **bold** indicate a sterile product

DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
<b>DOUBLE SWIVEL ELBOW</b>	-	<b>2401438-5</b>	<b>14 Fr/4.6mm</b>	<b>35cm</b>	<b>green</b>	<b>20/case</b>
Tracheostomy MDI Port						
<b>DOUBLE SWIVEL ELBOW</b>	-	<b>240168-5</b>	<b>16 Fr/5.3mm</b>	<b>58.5cm</b>	<b>orange</b>	<b>20/case</b>
Endotracheal MDI Port						
<b>DOUBLE SWIVEL ELBOW</b>	-	<b>2401638-5</b>	<b>16 Fr/5.3mm</b>	<b>35cm</b>	<b>orange</b>	<b>20/case</b>
Tracheostomy MDI Port						

All codes in **bold** indicate a sterile product

## BALLARD\* MULTI-ACCESS PORT CLOSED SUCTION SYSTEM – WET PAK\* DOUBLE SWIVEL ELBOWS

- BALLARD\* Technology
- Built-in Cleaning Chamber
- Rotational Manifold with Multiple Ports (Suctioning Port & Alternate Therapy for Bronchoscopy or BALLARD\* Mini-BAL sampling catheter)
- Tethered Sleeve
- Tethered Cap
- Available for Endotracheal and Tracheostomy Tubes
- Available with Metered Dose Inhaler (MDI) Port
- Contains in addition 12 Saline Vials
- Gamma Sterilised: 3 Years Shelf Life - Saline Vials: aseptically filled
- Suction Catheter: Validated for a maximum Usage Time of 72 Hours
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal	<b>FSQ1105</b>	<b>240126-5</b>	<b>12 Fr/4mm</b>	<b>58.5cm</b>	<b>white</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Tracheostomy	<b>FSQ1109</b>	<b>2401236-5</b>	<b>12 Fr/4mm</b>	<b>35cm</b>	<b>white</b>	<b>20/case</b>
<b>WET PAK DOUBLE SWIVEL ELBOW</b> Endotracheal MDI Port	<b>FSQ1113</b>	<b>2401286-5</b>	<b>12 Fr/4mm</b>	<b>58.5cm</b>	<b>white</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal	<b>FSQ1106</b>	<b>240146-5</b>	<b>14 Fr/4.6mm</b>	<b>58.5cm</b>	<b>white</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Tracheostomy	<b>FSQ1110</b>	<b>2401436-5</b>	<b>14 Fr/4.6mm</b>	<b>35cm</b>	<b>green</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal MDI Port	<b>FSQ1114</b>	<b>2401486-5</b>	<b>14 Fr/4.6mm</b>	<b>58.5cm</b>	<b>green</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Tracheostomy MDI Port	<b>FSQ1118</b>	<b>24014386-5</b>	<b>14 Fr/4.6mm</b>	<b>35cm</b>	<b>green</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal	<b>FSQ1107</b>	<b>240166-5</b>	<b>16 Fr/5.3mm</b>	<b>58.5cm</b>	<b>orange</b>	<b>20/case</b>


All codes in **bold** indicate a sterile product

DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
<b>WET PAK* DOUBLE SWIVEL ELBOW</b>	<b>FSQ1111</b>	<b>2401636-5</b>	<b>16 Fr/5.3mm</b>	<b>35cm</b>	<b>orange</b>	<b>20/case</b>
Tracheostomy						
<b>WET PAK* DOUBLE SWIVEL ELBOW</b>	<b>FSQ1115</b>	<b>2401686-5</b>	<b>16 Fr/5.3mm</b>	<b>58.5cm</b>	<b>orange</b>	<b>20/case</b>
Endotracheal MDI Port						
<b>WET PAK* DOUBLE SWIVEL ELBOW</b>	<b>FSQ1119</b>	<b>24016386-5</b>	<b>16 Fr/5.3mm</b>	<b>35cm</b>	<b>orange</b>	<b>20/case</b>
Tracheostomy MDI Port						

All codes in **bold** indicate a sterile product

## BALLARD\* MULTI-ACCESS PORT CLOSED SUCTION REPLACEMENT CATHETER

- BALLARD\* Technology
- Built-in Cleaning Chamber
- Gamma Sterilised: 3 Years Shelf Life
- Tethered Sleeve
- Tethered Cap
- Available for Endotracheal and Tracheostomy Tubes
- Suction Catheter: Validated for a maximum Usage Time of 72 Hours
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
<b>ENDOTRACHEAL</b> Replacement Catheter	<b>FSQ1121</b>	<b>240125-5</b>	<b>12 Fr/4mm</b>	<b>58.5cm</b>	<b>white</b>	<b>20/case</b>
<b>TRACHEOSTOMY</b> Replacement Catheter	<b>FSQ1125</b>	<b>2401253-5</b>	<b>12 Fr/4mm</b>	<b>35cm</b>	<b>white</b>	<b>20/case</b>
<b>ENDOTRACHEAL</b> Replacement Catheter	<b>FSQ1122</b>	<b>240145-5</b>	<b>14 Fr/4.6mm</b>	<b>58.5cm</b>	<b>green</b>	<b>20/case</b>
<b>TRACHEOSTOMY</b> Replacement Catheter	<b>FSQ1126</b>	<b>2401453-5</b>	<b>14 Fr/4.6mm</b>	<b>35cm</b>	<b>green</b>	<b>20/case</b>
<b>ENDOTRACHEAL</b> Replacement Catheter	<b>FSQ1123</b>	<b>240165-5</b>	<b>16 Fr/5.3mm</b>	<b>58.5cm</b>	<b>orange</b>	<b>20/case</b>


All codes in **bold** indicate a sterile product

# BALLARD\* CLOSED SUCTION SYSTEMS

**Indication:** Aspiration catheter for removal of tracheal secretions of intubated / tracheotomized adult patients respecting the closed ventilation circuit.

## BALLARD\* CLOSED SUCTION SYSTEM – T-PIECES

- BALLARD\* Technology
- Tactile Sleeve for Better Feel
- Gamma Sterilised: 5 Years Shelf Life
- Locking Thumb Port
- Validated for a maximum Usage Time of 24 Hours
- 22mm x 22mm Flex Tube Connector
- Available for Endotracheal and Tracheostomy tubes
- Available with Directional Tip
- Available with Dual Lumen Catheter
- Does not contain Natural Rubber Latex nor DEHP




DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
<b>CLEAR T-PIECE</b> Endotracheal	<b>FSQ1051</b>	<b>2205-5</b>	<b>14 Fr/4.6mm</b>	<b>54cm</b>	<b>green</b>	<b>20/case</b>
<b>CLEAR T-PIECE</b> Endotracheal Directional Tip	-	<b>22059-5</b>	<b>14 Fr/4.6mm</b>	<b>56cm</b>	<b>green</b>	<b>20/case</b>
<b>CLEAR T-PIECE</b> Tracheostomy 15mm x 22mm Swivel Adapter	<b>FSQ1052</b>	<b>220135-5</b>	<b>14 Fr/4.6mm</b>	<b>30.5cm</b>	<b>green</b>	<b>20/case</b>

All codes in **bold** indicate a sterile product

## BALLARD\* CLOSED SUCTION SYSTEM – DOUBLE SWIVEL ELBOWS

- BALLARD\* Technology
- Tactile Sleeve for Better Feel
- Numbered Catheter Allowing the Clinician to Limit the Distance of Catheter Advance
- Gamma Sterilised: 5 Years Shelf Life
- Available for Endotracheal and Tracheostomy tubes
- Available with Dual Lumen Catheter
- Available with Directional Tip
- Connector 15mm Female to 22mm Female
- Validated for a maximum Usage Time of 24 Hours
- Does not contain Natural Rubber Latex nor DEHP




DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal Directional Tip	-	<b>21609-5</b>	<b>12 Fr/4mm</b>	<b>56cm</b>	<b>white</b>	<b>20/case</b>
<b>DOUBLE SWIVEL ELBOW</b> Tracheostomy	<b>FSQ313</b>	<b>22103-5</b>	<b>14 Fr/4.6mm</b>	<b>30.5cm</b>	<b>green</b>	<b>20/case</b>
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal Dual Lumen	<b>FSQ268</b>	<b>22107-5</b>	<b>14 Fr/4.6mm</b>	<b>54cm</b>	<b>green</b>	<b>20/case</b>
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal MDI Port	<b>FSQ1058</b>	<b>22108-5</b>	<b>14 Fr/4.6mm</b>	<b>54cm</b>	<b>green</b>	<b>20/case</b>
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal Directional Tip	-	<b>22109-5</b>	<b>14 Fr/4.6mm</b>	<b>56cm</b>	<b>green</b>	<b>20/case</b>
<b>DOUBLE SWIVEL ELBOW</b> Endotracheal	<b>FSQ1060</b>	<b>2260-5</b>	<b>16 Fr/5.3mm</b>	<b>54cm</b>	<b>orange</b>	<b>20/case</b>

All codes in **bold** indicate a sterile product



## BALLARD\* CLOSED SUCTION SYSTEM – WET PAK\* DOUBLE SWIVEL ELBOWS

- BALLARD\* Technology
- Gamma Sterilised: 3 Years Shelf Life - Saline Vials: aseptically filled
- Tactile Sleeve for Better Feel
- Contains in addition 12 Saline Vials
- Validated for a maximum Usage Time of 24 Hours
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal	<b>FSQ419</b>	<b>21606-5</b>	<b>12 Fr/4mm</b>	<b>54cm</b>	<b>white</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Tracheostomy	<b>FSQ418</b>	<b>216036-5</b>	<b>12 Fr/4mm</b>	<b>30.5cm</b>	<b>white</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal	<b>FSQ421</b>	<b>22106-5</b>	<b>14 Fr/4.6mm</b>	<b>54cm</b>	<b>green</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Tracheostomy	<b>FSQ420</b>	<b>221036-5</b>	<b>14 Fr/4.6mm</b>	<b>30.5cm</b>	<b>green</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal Directional Tip	<b>FSQ321</b>	<b>221069-5</b>	<b>14 Fr/4.6mm</b>	<b>56cm</b>	<b>green</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Tracheostomy MDI Port	<b>FSQ1062</b>	<b>2210386-5</b>	<b>14 Fr/4.6mm</b>	<b>30.5cm</b>	<b>green</b>	<b>20/case</b>
<b>WET PAK* DOUBLE SWIVEL ELBOW</b> Endotracheal	<b>FSQ1063</b>	<b>22606-5</b>	<b>16 Fr/5.3mm</b>	<b>54cm</b>	<b>orange</b>	<b>20/case</b>

All codes in **bold** indicate a sterile product

# BALLARD\* CLOSED SUCTION SYSTEMS - ACCESSORIES

## BALLARD\* MULTI-ACCESS PORT CLOSED SUCTION SYSTEM - ACCESSORIES


- Gamma Sterilised




DESCRIPTION	NHSSC CODE	CODE	PACKAGING
BALLARD* SEAL CASSETTE	FSQ1128	24013-5	20/case

## BALLARD\* CLOSED SUCTION SYSTEM - ACCESSORIES

- Gamma Sterilised



118




DESCRIPTION	NHSSC CODE	CODE	PACKAGING
FLEX TUBE, 22MM X 22MM	-	<b>113</b>	50/case
Y-ADAPTER FOR SUCTION CANISTER	-	114	50/case
UNIT DOSE SALINE, 0.9% SODIUM CHLORIDE SOLUTION, PINK VIAL, 15ML, NONPYROGENIC, ASEPTICALLY FILLED	FSQ496	118	144/case
CONNECTOR 15MM FEMALE TO 22MM FEMALE	<b>FSQ613</b>	<b>1115</b>	50/case
SWIVEL ADAPTER, 15MM X 22MM, T-PIECE	-	<b>111</b>	50/case
FIXED ADAPTER, 15MM X 22MM, T-PIECE	<b>FSQ274</b>	<b>112</b>	50/case
METERED DOSE INHALER ADAPTER	-	<b>119</b>	20/case
BLUE ADAPTER FOR BALLARD* CLOSED SUCTION SYSTEM	-	50040	20/case

All codes in **bold** indicate a sterile product

## BALLARD\* CLOSED VENTILATOR CIRCUIT DRAIN

- Allows the Clearing of Ventilator Condensation Without Exposing the Caregiver or the Patient
- Vacuum Activated
- Maintains a Closed System
- Developed With the Safety of the Caregiver and Patient in Mind



DESCRIPTION	NHSSC CODE	CODE	PACKAGING
<b>BALLARD* CLOSED VENTILATOR CIRCUIT DRAIN ONLY</b>	-	131	<b>50/case</b>
<b>BALLARD* CLOSED VENTILATOR CIRCUIT DRAIN KIT:</b>	-	132	<b>25/case</b>
Circuit Drain 2,74 m Connector Tubing Y-Adapter			
<b>BALLARD* CLOSED VENTILATOR CIRCUIT DRAIN PROBE</b>	-	133	<b>50/case</b>
2,74 m Connector Tubing Y-Adapter			
<b>BALLARD* CLOSED VENTILATOR CIRCUIT DRAIN WITH 90 DEGREE CONNECTION</b>	-	134	<b>25/case</b>
2,74 m Connector Tubing Collapsible Connector Y-Adapter			

All codes in **bold** indicate a sterile product


# ENDOTRACHEAL TUBES - ADULTS

**Indication:** Airway management by nasal or oral intubation of the trachea in adult patients

- Oral/Nasal Magill
- Helps to prevent micro-aspiration
- Sealed airway at very low cuff pressure
- EO Sterilised: 5 Years Shelf Life
- Radio-opaque
- Anatomically-based Intubation Depth Mark
- Murphy Eye
- Does not contain Natural Rubber Latex

## MICROCUFF\* SUBGLOTTIC SUCTIONING ENDOTRACHEAL TUBES

- Available in sizes 7.0 to 9.0 mm (inner diameter)
- Cylindrical polyurethane ultra-thin cuff
- DEHP free



TUBE SIZE (INNER DIAMETER)	NHSSC CODE	CODE	PACKAGING
7.0 MM	<b>FD2117</b>	<b>13220</b>	1 case - 10 eaches/dispenser
7.5 MM	<b>FD2118</b>	<b>13221</b>	1 case - 10 eaches/dispenser
8.0 MM	<b>FD2119</b>	<b>13222</b>	1 case - 10 eaches/dispenser
8.5 MM	<b>FD2120</b>	<b>13223</b>	1 case - 10 eaches/dispenser
9.0 MM	<b>FD2121</b>	<b>13224</b>	1 case - 10 eaches/dispenser

## MICROCUFF\* ENDOTRACHEAL TUBES

- Available in sizes 5.0 to 10.0 mm (inner diameter)
- Cylindrical Cuff



TUBE SIZE (INNER DIAMETER)	NHSSC CODE	CODE	PACKAGING
5.0 MM	<b>FD1723</b>	<b>35210</b>	1 case - 10 eaches/dispenser
5.5 MM	<b>FD1724</b>	<b>35211</b>	1 case - 10 eaches/dispenser
6.0 MM	<b>FD1725</b>	<b>35212</b>	1 case - 10 eaches/dispenser
6.5 MM	<b>FD1726</b>	<b>35213</b>	1 case - 10 eaches/dispenser
7.0 MM	<b>FD1727</b>	<b>35214</b>	1 case - 10 eaches/dispenser
7.5 MM	<b>FD1728</b>	<b>35215</b>	1 case - 10 eaches/dispenser
8.0 MM	<b>FD1729</b>	<b>35216</b>	1 case - 10 eaches/dispenser
8.5 MM	<b>FD1730</b>	<b>35217</b>	1 case - 10 eaches/dispenser
9.0 MM	<b>FD1731</b>	<b>35218</b>	1 case - 10 eaches/dispenser
10.0 MM	-	<b>35220</b>	1 case - 10 eaches/dispenser

All codes in **bold** indicate a sterile product

# ORAL CARE SOLUTIONS - ADULTS

## BALLARD\* ORAL CARE KITS

**Indication:** Suction tools for removal of oral secretions of patients, that support the compliance to a comprehensive oral care protocol.



### BALLARD\* ORAL CARE ADVANCED KIT

- 2 Years Shelf Life

Kit includes:

- 1 Prep pack (incl. 1 covered Yankauer, 1 Y-adapter & 1 suction handle)
- 4 Suction Swab Packs (incl. 1 suction swab, 1 applicator swab & 1 empty solution cup)
- 2 Toothbrush Packs (incl. 1 suction toothbrush, 1 oral swab & 1 empty solution cup)

DESCRIPTION	NHSSC CODE	CODE	PACKAGING
BALLARD* ORAL CARE ADVANCED KIT	FSQ1071	25928	24/case

### BALLARD\* ORAL CARE MODERATE KIT

- 2 Years Shelf Life

Kit includes:

- 1 Prep Pack (incl. 1 covered Yankauer, 1 Y-adapter & 1 suction handle)
- 2 Suction Swab Packs (incl. 1 suction swab, 1 applicator swab & 1 empty solution cup)
- 1 Toothbrush Pack (incl. 1 suction toothbrush, 1 oral swab & 1 empty solution cup)

DESCRIPTION	NHSSC CODE	CODE	PACKAGING
BALLARD* ORAL CARE MODERATE KIT	FSQ1072	25930	32/case

All codes in **bold** indicate a sterile product

## BALLARD\* 24-HOUR ORAL CARE KIT (q4)

- 2 Years Shelf Life



Kit includes:

- 1 Prep Pack (incl. 1 covered Yankauer, 1 Y-adapter & 1 suction handle)
- 4 Suction Swab Procedure Packs (incl. 1 suction swab, 1 applicator swab, 1 empty solution cup)
- 2 Toothbrush Procedure Packs (incl. 1 suction toothbrush, 1 oral swab, 1 empty solution cup)
- 2 Catheter Packs (incl. 1 suction catheter)

DESCRIPTION	NHSSC CODE	CODE	PACKAGING
BALLARD* 24-HOUR ORAL CARE KIT q4	FSQ1070	99946	16/case

All codes in **bold** indicate a sterile product

## BALLARD\* MINI ORAL CARE KIT

- 2 Years Shelf Life
- For use with a 4-hour oral care protocol



Kit includes:

- 1 Prep Pack (reusable for 24H - incl. self-cleaning covered Yankauer)
- 2 Toothbrush Packs (incl. 1 suction toothbrush, 1 mini applicator swab, 1 empty solution cup & 1 mini suction catheter)
- 4 Suction Swab Packs (incl. 1 mini suction swab, 1 mini applicator swab & 1 empty solution cup)

DESCRIPTION	NHSSC CODE	CODE	PACKAGING
BALLARD* MINI ORAL CARE KIT	FDB953	970008-5	24/case


All codes in **bold** indicate a sterile product

# BALLARD\* ORAL CARE SUCTION KITS - COMPONENTS & SUCTION KITS

**Indication:** Suction tools for removal of oral secretions of patients, that support the compliance to a comprehensive oral care protocol.

## BALLARD\* ORAL CARE SUCTION KITS

- 2 Years Shelf Life




DESCRIPTION	NHSSC CODE	CODE	PACKAGING
<b>BALLARD* CONNECTION TUBING WITH SUCTION VALVE:</b>	FSQ718	122	<b>50/case</b>
1 suction valve, 2.74 m funnel tubing 1 suction line clamp, 1 Y-adaptor			

All codes in **bold** indicate a sterile product



## BALLARD\* ORAL CARE SUCTION COMPONENTS

- 5 Years Shelf Life (except 99785<sup>◊</sup>, 99786<sup>†</sup>, 99789<sup>‡</sup>: 2 Years Shelf Life)




DESCRIPTION	NHSSC CODE	CODE	PACKAGING
<b>BALLARD* READY CARE* ORAL CARE SUCTION PROBES</b>	FSQ719	123	<b>50/case</b>
<b>BALLARD* ORAL CARE SUCTION CATHETER FOR ADULTS</b>	FSQ373	12601	<b>25/case</b>
<b>BALLARD* ORAL CARE SUCTION TOOTHBRUSH</b>	FSQ374	12602	<b>25/case</b>
<b>BALLARD* ORAL CARE SUCTION SWAB, ANGLED-TIP</b>	-	12603	<b>25/case</b>
<b>BALLARD* ORAL CARE SUCTION TIP YANKAUER</b>	FSQ376	12604	<b>25/case</b>
<b>BALLARD* SUCTION VALVE AND Y-ADAPTER</b>	FSQ377	12606	<b>25/case</b>
<b>BALLARD* ORAL CARE SUCTION SWAB, WITH DENTIFRICE (MAY BE USED WITHOUT BALLARD* SUCTION VALVE)</b>	FSQ1888	12607	<b>25/case</b>
<b>BALLARD* SELF-CLEANING COVERED YANKAUER</b>	FSQ2400	99785 <sup>◊</sup>	<b>25/case</b>
<b>BALLARD* SUCTION HANDLE</b>	FSQ1284	99786 <sup>†</sup>	<b>25/case</b>
<b>BALLARD* ORAL CARE PREP PACK:</b>	FSQ1285	99789 <sup>‡</sup>	<b>40/case</b>
1 self-cleaning, covered Yankauer, 1 suction handle, Y-adapter			

All codes in **bold** indicate a sterile product

## BALLARD\* ORAL CARE SWABS

- Shelf Life for Swabs with Dentifrice: 2 years
- Shelf Life for Swabs without Dentifrice: 5 years



DESCRIPTION	NHSSC CODE	CODE	ITEMS/POUCH	POUCHES/CASE
<b>SINGLE WRAPPED</b>				
BALLARD* ORAL CARE SWAB	FSQ1881	12240	250	4
BALLARD* ORAL CARE SWAB, WITH DENTIFRICE	FSQ1884	12245	250	4
BALLARD* ORAL CARE SWAB, WITH DENTIFRICE	-	12246	250	2
BALLARD* ORAL CARE SWAB, PLASTIC HANDLE	FSQ1886	12254	250	4
BALLARD* ORAL CARE SWAB, WITH DENTIFRICE, PLASTIC HANDLE	FSQ1887	12255	250	4
<b>BULK PACKAGED</b>				
BALLARD* ORAL CARE SWAB	FSQ1882	12241	500	2
BALLARD* ORAL CARE SWAB	FSQ1883	12243	20	50
BALLARD* ORAL CARE SWAB, WITH DENTIFRICE	FSQ1885	12249	20	50
BALLARD* ORAL CARE SWAB, PLASTIC HANDLE	-	12251	20	50
BALLARD* ORAL CARE SWAB, WITH DENTIFRICE, PLASTIC HANDLE	-	12257	20	50

## BALLARD\* MINI ORAL CARE SWABS

- Shelf Life for Swabs: 5 years
- Individually Wrapped with 12 cm Plastic Handle

DESCRIPTION	NHSSC CODE	CODE	ITEMS/POUCH	POUCHES/CASE
BALLARD* MINI ORAL CARE SWABS	FDB952	122548	250	4


All codes in **bold** indicate a sterile product

# DIAGNOSTIC CATHETERS - ADULTS

## BALLARD\* MINI-BAL SAMPLING CATHETER

**Indication:** This product is intended for early, accurate diagnosis of VAP.

- Simple, effective non-bronchoscopic bronchoalveolar lavage at the bedside
- EO Sterilised: 3 Years Shelf Life
- Directional Tip for R/L Insertion
- Includes an access port elbow, adapter for seal cassette, seal cassette, stepped connector, oxygen line
- Atraumatic "Mushroom" Tip
- Does not contain Natural Rubber Latex



DESCRIPTION	NHSSC CODE	CODE	PACKAGING
BALLARD* MINI-BAL SAMPLING CATHETER - 16 FR OUTER/12 FR INNER TELESCOPING CATHETER	<b>FSQ1914</b>	<b>142</b>	5/case
BALLARD* MINI-BAL SAMPLING CATHETER - 13 FR OUTER/8 FR INNER TELESCOPING CATHETER	<b>FSQ1913</b>	<b>143</b>	5/case

All codes in **bold** indicate a sterile product

# CLOSED SUCTION SYSTEMS - NEONATES/PAEDIATRICS

**Indication:** Aspiration catheter for removal of tracheal secretions of intubated neonatal/paediatric patients respecting the closed ventilation circuit.

## BALLARD\* CLOSED SUCTION SYSTEM FOR NEONATES/PAEDIATRICS - ELBOWS

- BALLARD\* Technology
- Gamma Sterilised: 5 Years Shelf Life
- Locking Thumb Port
- Contains Elbow Adapters
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
ELBOW ENDOTRACHEAL TUBE ADAPTERS (2.5mm, 3mm, 3.5mm)	<b>FSQ499</b>	<b>206-5</b>	<b>6 Fr/2mm</b>	<b>30.5cm</b>	<b>light green</b>	<b>20/case</b>
ELBOW ENDOTRACHEAL TUBE ADAPTERS (2.5mm, 3mm, 3.5mm, 4mm)	<b>FSQ1065</b>	<b>207-5</b>	<b>7 Fr/2.3mm</b>	<b>35.5cm</b>	<b>ivory</b>	<b>20/case</b>
ELBOW ENDOTRACHEAL TUBE ADAPTERS (3mm, 3.5mm, 4mm)	<b>FSQ311</b>	<b>208-5</b>	<b>8 Fr/2.6mm</b>	<b>35.5cm</b>	<b>light blue</b>	<b>20/case</b>
ELBOW ENDOTRACHEAL TUBE ADAPTERS (4.5mm, 5mm, 5.5mm)	<b>FSQ498</b>	<b>210-5</b>	<b>10 Fr/3.3mm</b>	<b>40.5cm</b>	<b>black</b>	<b>20/case</b>

All codes in **bold** indicate a sterile product

## BALLARD\* CLOSED SUCTION SYSTEM FOR NEONATES/PAEDIATRICS - "Y"

- BALLARD\* Technology
- Gamma Sterilised: 5 Years Shelf Life
- Locking Thumb Port
- Contains Elbow and Y-adapters
- Does not contain Natural Rubber Latex nor DEHP



DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
"Y" ENDOTRACHEAL TUBE ADAPTERS (2mm, 2.5mm)	FSQ300	195-5	5 Fr/1.6mm	30.5cm	grey	20/case
"Y" ENDOTRACHEAL TUBE ADAPTERS (2.5mm, 3mm, 3.5mm)	FSQ301	196-5	6 Fr/2mm	30.5cm	light green	20/case
"Y" ENDOTRACHEAL TUBE ADAPTERS (3mm, 3.5mm, 4mm)	FSQ1066	197-5	7 Fr/2.3mm	30.5cm	ivory	20/case
"Y" ENDOTRACHEAL TUBE ADAPTERS (3mm, 3.5mm, 4mm)	FSQ308	198-5	8 Fr/2.6mm	30.5cm	light blue	20/case
"Y" ENDOTRACHEAL TUBE ADAPTERS (4.5mm, 5mm, 5.5mm)	FSQ290	1910-5	10 Fr/3.3mm	40.5cm	black	20/case
"Y" ENDOTRACHEAL TUBE ADAPTERS (5.5mm, 6mm, 6.5mm)	FSQ291	1912-5	12 Fr/4mm	40.5cm	white	20/case

All codes in **bold** indicate a sterile product

## BALLARD\* MULTI-ACCESS CATHETER FOR NEONATES/PAEDIATRICS


**Indication:** A multi-access catheter designed to deliver surfactants to the distal airway

- BALLARD\* Technology
- Gamma Sterilised: 5 Years Shelf Life
- Contains Closure Cap and Y-adapter
- Does not contain Natural Rubber Latex
- The Catheter is DEHP free



DESCRIPTION	NHSSC CODE	CODE	OUTER DIAMETER	LENGTH	ISO COLOUR	PACKAGING
MAC CATHETER ADAPTER (2.5mm, 3mm, 3.5mm)	FSQ403	1900	5 Fr/1.6mm	30.5cm	white	5/case
MAC CATHETER ADAPTER (2mm)	FSQ1067	1920	5 Fr/1.6mm	30.5cm	white	5/case
MAC CATHETER ADAPTER (4mm)	FSQ1068	1940	5 Fr/1.6mm	30.5cm	white	5/case
TRANSPORT KIT ("Y" Plug and Cap)	FSQ1069	19000				50/case

## BALLARD\* CLOSED SUCTION SYSTEM FOR NEONATES/PAEDIATRICS - ACCESSORIES



DESCRIPTION	NHSSC CODE	CODE	PACKAGING
ENDOTRACHEAL TUBE ADAPTER 2MM WITH LAVAGE PORT	-	2-020	20/case
ENDOTRACHEAL TUBE ADAPTER 2.5MM WITH LAVAGE PORT	-	2-025	20/case
ENDOTRACHEAL TUBE ADAPTER 3MM WITH LAVAGE PORT	-	2-030	20/case
ENDOTRACHEAL TUBE ADAPTER 3.5MM WITH LAVAGE PORT	-	2-035	20/case
ENDOTRACHEAL TUBE ADAPTER 4MM WITH LAVAGE PORT	-	2-040	20/case
ENDOTRACHEAL TUBE ADAPTER 4.5MM WITH LAVAGE PORT	-	2-045	20/case
MANUAL PERCUSSOR FOR NEONATAL PULMONARY USE	-	11000	12/case


All codes in **bold** indicate a sterile product

# ENDOTRACHEAL TUBES - PAEDIATRICS

**Indication:** Airway management by nasal or oral intubation of the trachea in paediatric patients

- Available in Magill & RAE / preformed types
- Available in sizes 3.0 to 7.0 mm (inner diameter)
- EO Sterilised: 5 Years Shelf Life
- Designed for ideal anatomic cuff placement
- Sealed airway at ultra low cuff pressure
- Does not contain Natural Rubber Latex


## MICROCUFF\* ENDOTRACHEAL TUBE FOR PAEDIATRIC, ORAL/NASAL MAGILL



TUBE SIZE (INNER DIAMETER)	NHSSC CODE	CODE	AGE/WEIGHT YEARS/KG	PACKAGING
3.0 MM	FDG1209	<b>35111</b>	Term ≥ 3kg to < 8 months	10/dispenser
3.5 MM	FDG1210	<b>35112</b>	8 months < 2 years	10/dispenser
4.0 MM	FDG1211	<b>35113</b>	2 to < 4 years	10/dispenser
4.5 MM	FDG1212	<b>35114</b>	4 to < 6 years	10/dispenser
5.0 MM	FDG1213	<b>35115</b>	6 to < 8 years	10/dispenser
5.5 MM	FDG1214	<b>35116</b>	8 to < 10 years	10/dispenser
6.0 MM	FDG1215	<b>35117</b>	10 to < 12 years	10/dispenser
6.5 MM	FDG1216	<b>35118</b>	12 to < 14 years	10/dispenser
7.0 MM	FDG1217	<b>35119</b>	14 to < 16 years	10/dispenser

All codes in **bold** indicate a sterile product

## MICROCUFF\* ENDOTRACHEAL TUBE FOR PAEDIATRIC, ORAL CURVED



TUBE SIZE (INNER DIAMETER)	NHSSC CODE	CODE	AGE/WEIGHT YEARS/KG	PACKAGING
3.0 MM	FDG1577	35161	Term $\geq$ 3kg to < 8 months	10/dispenser
3.5 MM	FDG1575	35162	8 months < 2 years	10/dispenser
4.0 MM	FDG1569	35163	2 to < 4 years	10/dispenser
4.5 MM	FDG1570	35164	4 to < 6 years	10/dispenser
5.0 MM	FDG1571	35165	6 to < 8 years	10/dispenser
5.5 MM	FDG1576	35166	8 to < 10 years	10/dispenser
6.0 MM	FDG1572	35167	10 to < 12 years	10/dispenser
6.5 MM	FDG1573	35168	12 to < 14 years	10/dispenser
7.0 MM	FDG1574	35169	14 to < 16 years	10/dispenser

All codes in **bold** indicate a sterile product



# ORAL CARE SOLUTIONS - NEONATES/PAEDIATRICS

## BALLARD\* ORAL CARE SUCTION CATHETERS FOR NEONATES/PAEDIATRICS

- Does not contain Natural Rubber Latex
- Packaging: 5 pouches of 50 units

DESCRIPTION	NHSC CODE	CODE	PACKAGING
<b>5 BALLARD* ORAL CARE SUCTION CATHETERS FOR NEONATES/PAEDIATRICS</b> (6 Fr - 2 mm)	-	1221	<b>50/case</b>
<b>5 BALLARD* ORAL CARE SUCTION CATHETERS FOR NEONATES/PAEDIATRICS</b> (8 Fr - 2.6 mm)	-	1222	<b>50/case</b>
<b>5 BALLARD* ORAL CARE SUCTION CATHETERS FOR NEONATES/PAEDIATRICS</b> (10 Fr - 3.3 mm)	FSQ546	1223	<b>50/case</b>

All codes in **bold** indicate a sterile product

# INDEX BY CODE

CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE
2-020	30	12243	26	35118	31	240165-5	14
2-025	30	12245	26	35119	31	240166-5	12
2-030	30	12246	26	35161	32	240168-5	11
2-035	30	12249	26	35162	32	970008-5	23
2-040	30	12251	26	35163	32	2210386-5	17
2-045	30	12254	26	35164	32	2270135	8
111	18	12255	26	35165	32	2270156-5	9
112	18	12257	26	35166	32	2271018-5	5
113	18	12601	25	35167	32	2271053-5	8
114	18	12602	25	35168	32	2271218-5	5
118	18	12603	25	35169	32	2271418-5	6
119	18	12604	25	35210	20	2271516-5	9
122	24	12606	25	35211	20	2271603-5	5
123	25	12607	25	35212	20	2271613-5	6
131	19	13220	20	35213	20	2271616-5	7
132	19	13221	20	35214	20	2271618-5	6
133	19	13222	20	35215	20	2401236-5	12
134	19	13223	20	35216	20	2401253-5	14
142	27	13224	20	35217	20	2401286-5	12
143	27	19000	30	35218	20	2401436-5	12
195-5	29	21606-5	17	35220	20	2401438-5	11
196-5	29	21609-5	16	50040	18	2401453-5	14
197-5	29	22059-5	15	99785*	25	2401486-5	12
198-5	29	22103-5	16	99786*	25	2401636-5	13
206-5	28	22106-5	17	99789*	25	2401638-5	11
207-5	28	22107-5	16	99946	22	2401686-5	13
208-5	28	22108-5	16	122548	26	22012356-5	9
210-5	28	22109-5	16	216036-5	17	22701356-5	9
227-5	5	22606-5	17	220135-5	15	22712183	5
1115	18	22703-5	6	221036-5	17	22712186-5	7
1221	33	22715-5	8	221069-5	17	22714183-5	6
1222	33	22716-5	5	227036-5	7	22714186-5	7
1223	33	24010-5	10	227101-5	5	22716036-5	7
1900	30	24013-5	18	227105-5*	8	22716183-5	6
1910-5	29	24014-5	10	227135-5	8	22716186-5	7
1912-5	29	25928	21	227161-5	6	24001238-5	10
1920	30	25930	21	227166-5	7	24014386-5	12
1940	30	35111	31	240121-5	10	24016386-5	13
2205-5	15	35112	31	240125-5	14	227121863-5	7
2260-5	16	35113	31	240126-5	12	227141863-5	7
2276-5	7	35114	31	240128-5	10	227161863-5	7
11000	30	35115	31	240145-5	14		
12240	26	35116	31	240146-5	12		
12241	26	35117	31	240148-5	10		



---

**AVANOS**

From the UK, please call Avanos  
at 0800 917 65 85 - fax 0800 169 02 35  
or send an email to [customerservice.uk.ie@avanos.com](mailto:customerservice.uk.ie@avanos.com)

From all other countries, please call Avanos  
at +32 2 700 68 51 - fax +32 2 711 26 91  
or send an email to [customerservice.export@avanos.com](mailto:customerservice.export@avanos.com)

[www.avanos.co.uk](http://www.avanos.co.uk)