



## AccuPAP™

### Lung Expansion Therapy in the Palm of Your Hand

- Provides positive airway pressure for the prevention and treatment of atelectasis.
- Optimizes delivery of aerosolized medications.
- PEEP independent of flow.
- Four PEEP Settings: 5/10/15/20cmH<sub>2</sub>O.

## AccuPAP™- Positive Airway Pressure Device,

Part No.	Product Description	Case/Qty.
313-6000	AccuPAP with Mouthpiece and 7' Tubing	10
313-6001	AccuPAP with Mouthpiece, 7' Tubing, and DPM™ (Disposable Pressure Manometer)	10



## DPM™- Disposable Pressure Manometer

Part No.	Product Description	Case/Qty.
313-5558	DPM (Disposable Pressure Manometer)	10

