

WatchBP Cuffs

The cuff is as important as the blood pressure monitor itself. Using the right size cuff helps to prevent erroneous blood pressure measurements.




WatchBP Cuffs

- ✓ Latex-free and PVC-free.
- ✓ Washable cuffs available in both nylon and cotton.
- ✓ Disposable cuffs available (for single patient use).

WatchBP Disposable cuffs

Upper arm – Single patient use (For WatchBP Office series)


S




M




L



Tube with cuff connector




50 cm tube with cuff connector




100 cm tube with cuff connector

WatchBP Office Cuffs

Upper arm – Soft




Cuff connector




For Office AFIB, Office ABI

Upper arm – Rigid

M-L
 With air tube 130 cm for Office AFIB, Office ABI




Cuff connector




For Office AFIB, Office ABI

Ankle – Soft



Tube with cuff connector



170 cm tube with cuff connector

