

**PHILIPS**

RESPIRONICS

***Wisp***

Family nasal mask

Rediscover  
**dreams**



## Wisp family nasal masks serve infants to adults for continuity and consistency in their therapy experience.

A range of cushion and headgear sizes are available to accommodate up to 98% of your patients.

Wisp pediatric is designed for infants weighing 5 kg to 18 kg (approx. 7 years old).

Wisp youth contains three cushion sizes (Petite, Small/Medium and Large) and fits 97% of children from 7 to 18 years old. Wisp youth is intended for children weighing 18 kg to 30 kg.

Wisp adult also contains three cushion sizes: Small/Medium, Large and Extra Large, and will fit 98% of adults. Wisp adult is intended for adults weighing 30 kg and up.



Petite

Small/Medium

Large

Extra Large

# Wisp pediatric

## Unique design for unique patient needs

Wisp pediatric cushion, frame and headgear are designed to fit the specific sizing and bone structure of your pediatric patients.

- Open line of sight allows kids to comfortably wear the mask during their normal daily activities.
- Cushion design features Auto Seal technology, delivering minimal contact and a snug seal.
- 360° rotating elbow
- Click fit headgear clips, with adjustable straps for quick and easy adjustments.



# Wisp adult/youth

## 'Tip-of-the-nose' nasal mask

Wisp's clever Auto Seal technology accommodates your patient's every move, for a superior seal and a great night's sleep.


- Headgear with elastic crown strap and integrated tube holder. Self-adjusting for a more natural fit
- 360° rotating elbow
- Contemporary frame options to fit your patient's style: Choose from a silicone/clear frame and a reversible fabric frame - soft suede or silky sateen.



**Magnetic headgear clips**



**Wisp silicone frame**



## Leak

### Correction Dial

for less intrusive mask adjustment. If leaks occur, turn the dial left and right, gently nudging the cushion closer to the face

### Optional tube wrap

a tubing management loop enables the tube to be positioned over the head and away from the front of the body



## 'Tip-of-the-nose'

**cushion design, featuring Auto Seal technology**

Delivers a minimal contact experience with an exceptional seal

### Magnetic headgear clips\* with adjustable straps

enhance ease of use for patients with limited dexterity

\*Intended for patients weighting 66 lbs / 30 kgs and up. To be used with Wisp adult nasal mask only.



## Supporting you and your patients

There are a range of materials available, including:

- Adult/youth and pediatric fitting guides
- Instructional videos and 360° view on the product web page
- Supporting tools for parents and children, such as the Wisp pediatric Caregiver guide, Storybook and the Jacky Giraffe adventure cartoon.

**For more information, please visit**  
[www.philips.com/respironics](http://www.philips.com/respironics)

<b>Part no.</b>	<b>Description</b>
1104961	Wisp pediatric with exhalation port and child size headgear
1109308	Wisp youth with exhalation port, fabric frame and small size headgear
1094060	Wisp adult with exhalation port, silicone frame and standard size headgear
1094061	Wisp adult with exhalation port, fabric frame and standard size headgear
1094062	Wisp adult without exhalation port, silicone frame and standard size headgear



SJ 9/20/17 PN499439  
[www.philips.com/respironics](http://www.philips.com/respironics)

© 2017 Koninklijke Philips N.V. All rights reserved.  
Specifications are subject to change without notice.