


# Supporting Your Telehealth Needs

Telehealth solutions are becoming increasingly critical in the global healthcare environment. RescueNet Live from ZOLL allows clinicians to remotely view real-time data from the X Series® Advanced monitor/defibrillator screen via a cellular or WiFi connection.



RescueNet Live allows providers access to:

- Top 4 waveforms\*
- Heart rate
- Non-invasive blood pressure
- EtCO<sub>2</sub> values
- Masimo rainbow® SET® parameters
- Invasive blood pressure
- Temperature

Providers will also have visibility into:

- Trends
- 12-lead snapshots
- Alarms
- Patient demographics
- Gallery view of multiple streams
- Multiple users connecting to a specific stream

RescueNet Live offers effective decision support to help EMS organizations meet their telehealth needs.

As telehealth programs continue to help improve patient care, providers can rely on RescueNet Live to share the real-time feedback and data necessary to provide optimal care in every situation.

For more information, visit [zoll.com/rescuenetlive](https://zoll.com/rescuenetlive)

\*Available in version 2.36 and later.  
2.36 is not currently available in North America.