

O2-MAX™ Trio

Emergency CPAP with the Ability to Increase your FiO2

- The O2-MAX uses a 50 PSI port, and produces flows of up to 140 LPM.
- It has an approximate 30%, 60%, or 90% FiO2.
- It includes a 3-SET™ option for PEEP of 5, 7.5, and 10 cm H2O, and a 5-SET™ option of 5, 7.5, 10, 12.5, and 15 cm H2O.
- It offers integrated nebulization without interrupting the delivered flow to the patient.



O2-MAX™ Trio system with BiTrac™ ED Mask and Expandable Anti-Asphyxia Circuit - *no nebulizer canister included*

Part No.	PEEP Setting	Size	Case/Qty.
313-8007X	3-SET™ PEEP	Adult Medium	10
313-8017X	3-SET™ PEEP	Adult Large	10
313-8008X	5-SET™ PEEP	Adult Medium	10
313-8018X	5-SET™ PEEP	Adult Large	10
313-8003X	10cm Fixed PEEP	Adult Medium	10

O2-MAX™ Trio system with BiTrac™ ED Mask and Expandable Anti-Asphyxia Circuit - *nebulizer canister and Quik-Connect™ included*

Part No.	PEEP Setting	Size	Quik-Connect™	Case/Qty.
313-8007XN-1	3-SET™ PEEP	Adult Medium	Ohmeda	10
313-8017XN-1	3-SET™ PEEP	Adult Large	Ohmeda	10
313-8007XN-2	3-SET™ PEEP	Adult Medium	Chemtron	10
313-8007XN-13	3-SET™ PEEP	Adult Medium	AFNOR & DIN	10
313-8017XN-13	3-SET™ PEEP	Adult Large	AFNOR & DIN	10

Individual Headgear with clips

Part No.	Size	Case/Qty.
313-7028	Adult Small	10
313-7030	Adult Large	10

