





## Control-Cric™

### Adding Control to an Otherwise Chaotic Procedure

- Surgical cricothyrotomy device for use in an emergency airway setting.
- Used when airway obstruction has occurred, or when a patient cannot be intubated or ventilated.
- Reliable and intuitive device, with a three step compact solution.
- Allows the clinician to maintain placement with an integrated sliding tracheal hook that sits within the knife.
- Provides tactile feedback of the tracheal rings, before placement of the tube by using an integrated stylet.



### Control-Cric™ System- Sterile kit includes a Cric-Knife™, Cric-Key™, Neck Strap, Syringe, and extension tubing.

Part No.	Product Description	Case/Qty.	
313-4567	Control-Cric System (Sterile)	5	
			

### Control-Cric™ Training System Non-Sterile kit includes a Cric-Knife™, Cric-Key™, Neck Strap, Syringe, and extension tubing. \*Not For Human Use

Part No.	Product Description	Case/Qty.	
313-4568	Training Control-Cric System (Non-Sterile)	10	
